

## PORTABLE RAMAN SPECTROMETER PART No. HHL-R550

- Used in the fields of medicine, environmental protection, food, public security, etc.
- Standard database, supports users to add and modify spectrograms
- Analysis of mixture algorithmNon-destructive and rapid detection
- Detection records can be exported
- Integrated lithium battery, which can work continuously for 4-6 hours when fully charged



HHL-R550

## SPECIFICATION

SPECIFICATION		
Exciting Laser	central wavelength	785±1nm
	FWHM	0.08nm
	output power	0~500mW
	power stability	σ/μ<±0.2%
Detector	type	ultra-high sensitivity, quick cooling CCD
	refrigeration temperature	14°F
	detecting range	200~1100nm
	effective pixels	2048×64
	dynamic range	50000:1
	pixel size	14μm×14μm
	full well capacity	300Ke <sup>-</sup>
	sensitivity	QE>40%, 6.5μV/e <sup>-</sup>
Optical parameters	spectral range	250~2700cm <sup>-1</sup>
	resolution	6cm <sup>-1</sup>
	entrance slit	50μm
	focusing	98mm
Raman probe	operating distance	6mm
	rayleigh scattering resistance	OD>8
	numerical aperture	0.3
	aperture	7mm
Spectral stability		σ/μ<0.5% (8 hours)
Temp stability		spectral shift<1cm <sup>-1</sup> (50~104°F)
Intensity repeatability		≤5%
Frequency shift repeatability		≤1cm <sup>-1</sup>
Frequency shift accuracy		≤1cm <sup>-1</sup>
Integration time		4ms~120s
Operating system		Android
Interface		WIFI, USB
Working environment		41~104°F, ≤95%RH, no condensation
Storage environment		-4~113°F, ≤95%RH, no condensation
Power supply		19V rechargeable lithium battery
Dimension (L×W×H)		16.93x15.35x6.69"
Weight		20.50 <b>l</b> b

## STANDARD DELIVERY

Main unit	1 pc
Liquid sampling accessory	1 pc
Solid sampling accessory	1 pc
Power adapter	1 pc

## **OPTIONAL ACCESSORY**

Test adjustment frame	HHL-R550-TS
_ , . , . , . ,	
Test adjustment frame	HHL-R550-TS

applications



pharmaceuticals





