

# AUTOMATIC LIQUID DROPPER



MLP-PLD101-U



MLP-PLD101S-U

- Adopting advanced micro-controller, compact structure, reliable operation, strong anti-interference ability
- Suitable for various types of grinding and polishing machines, the polishing liquid is automatically dripped between the abrasive disk and the sample tray to achieve the purpose of polishing
- Reduced polishing time, increased polishing efficiency



Polishing fluid (optional)

## SPECIFICATION

Code	MLP-PLD101-U (single type)	MLP-PLD101S-U (combination type)
Dropping rate	0.3~10ml/min	
Number of steps	10-speed	
Tank volume	750ml	
Power supply	AC 110V, 60HZ	
Dimension(W×H×D)	5.67×9.09×7.76"	22.84×9.25×8.35"
Weight	6.51lb	42.99lb

## STANDARD DELIVERY

Main unit	1 pc
Diamond polishing liquid	1 pc (W2.5, 500ml)

## OPTIONAL DELIVERY

Diamond polishing liquid	MLP-DLA□□*, MLP-DLB□□*, MLP-DLC□□*, MLP-DLD□□*, MLP-DLE□□*, MLP-DLF□□*
Aluminum oxide polishing liquid	MLP-LXA□□*, MLP-LXB□□*
Silicon dioxide polishing liquid	MLP-GF1, MLP-GF2, MLP-GF3
Polishing Lubricant	MLP-RH1, MLP-RH2, MLP-RH3

\*□□ is granularity specification, for example, code MLP-DLA40 stands for 40um